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13~15

Introduction

Thank you for purchasing the Naipo Folding Portable Massage Chair. Please read this user manual carefully and keep it in a safe place for future reference. If you need any assistance, please contact our support team with your product model number and Amazon order number.

Included

1 x Portable Massage Chair
1 x User Manual
1 x Carrying Bag
1 x Headrest
1 x Seat cushion

Product Specs

Model	MGCHR-001H
Name	Naipo Folding Portable Massage Chair

Important Safeguards

1. Maintain the chair regularly. Make sure all screws are secure and that the components are complete without damage.

2. Do not overload the chair beyond the recommended weight limit-the combined weight of the client and the maximum amount of downward pressure being applied by the therapist.

3. Do not clean the product surface with alcohol-based cleaner, otherwise it will cause damage to the product; Clean with soapy water and then wipe dry.

4. When using the chair, please make sure that all four legs are placed stably on a flat, solid surface.

Installation Guide

Unpacking
Take out the massage chair from the package and keep upright.

Adjusting
1. Grab the front frame and pull it outward, and let the knee pads fall down on the front frame. After the massage chair is unfolded, install the seat cushion, then adjust it to an appropriate angle and height, with seat support tube being placed securely on the corresponding position of the support slot.

User tip: the seat cushion is included separately in the package and needs to be manually installed, and it comes with the quick release device. Accessories like screws and quick release lever are included on the conjunction area where the seat support plate is installed.

2. The headrest is also included separately and needs to be installed manually, insert the headrest tube inside the headrest support tube on the massage chair, with the lightness being adjustable through the knob. To move the headrest up or down, please loosen the knob on the chest pad support plate and then pull (or push) the headrest tube. To adjust the headrest angle, please unlock the quick release device, then lock again by flipping the lever reversely after the adjustment.

- 01 -

Contents

1. Move the chest pad and headrest to the lowest position.
2. Loosen the armrest and let it drop to the most vertical position.
3. Fold the headrest forward fully.
4. Lift the seat and pull up the knee pads until the rear frame moves close to the front frame.
5. Store in the carrying bag.

Warranty & Customer Support

If you need help regarding product issues, after-sales service, or other relevant problems, please contact us via the email on the cover.

*Note: Products purchased directly from Naipo will enjoy our after-sales service. If the purchases are made from other selling teams, we recommend you contact them for warranty issues.

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